



EXOFORM® KNEE IMMOBILIZER

The Exoform Knee Immobilizer is the most adjustable knee immobilizer available. Össur's unique Slide-to-Size™ approach enables the care provider to fit the widest variety of patients. This allows you to reduce your inventory. The Exoform's unique posterior stay and dual-cuff design delivers greater immobilization compared to the old-style softgood approach.

- Slide-to-size adjustment enables you to provide the right fit for the largest variety of patient lengths and widths
- Dual-cuff design leads to the greatest immobilization
- Popliteal pad provides superior fit and limits migration
- Perforated padded liner for breathability and comfort
- Cool foam version can be worn directly against the skin for greater comfort on warm days
- Full foam version provides warmth and compression
- Uncompromising durability
- Easy to fit patient quickly

PART#	FOAM	LEG CIRC.	LENGTH	U/M
222100	Cool	up to 28"	18-24"	1 ea
222000	Full	up to 28"	18-24"	1 ea
222050	XXXL-Full	up to 38"	18-24"	1 ea

Replacement Softgoods

PART#	FOAM
222101	Cool
222001	Full
222051	XXXL-Full

U.S. and Foreign Patents Pending.

UNIVERSAL 3-PANEL KNEE IMMOBILIZER

- Two-sided brushed nylon, laminated to breathable foam
- Three piece design for precise fit
- Elastic patellar strap
- Velcro straps with loop closures
- Lightweight rigid posterior stays with adjustable medial and lateral stays.

PART#	LENGTH
200010	14"
200016	16"
200118	18"
200020	20"
200022	22"
200030	24"

Add the healing properties of safe Hot/Cold Therapy to your Knee Immobilizer.

Use part# GB7

