

TRAINER OTS

This cost-effective, single-upright brace is the perfect option to treat moderate MCL injuries, or mild combined ligament instabilities. With a diagonal Dynamic Force Strap¹ that tightens when the knee is extended, Trainer resists valgus forces and minimizes axial rotation. Its carbon steel range-of-motion hinge² and injection-molded carbon shells³ are extremely durable, and an Adjustable Tibial Strap⁴ provides tibial capture and improved suspension. The Trainer is an excellent option for use in activities where a medial hinge may impede performance, or a padded shell covering is required.



INDICATIONS

- Mild to moderate ligament laxity and acute management of Grade I and Grade II MCL sprains
- Recommended for low to medium contact/impact

PART#	LEG	SIZE	CIRC. ¹
B-237500001	Right	XSmall	13.5-15.5"
B-237600001	Left	XSmall	13.5-15.5"
B-237500002	Right	Small	15.5-17.5"
B-237600002	Left	Small	15.5-17.5"
B-237500003	Right	Medium	17.5-20.5"
B-237600003	Left	Medium	17.5-20.5"
B-237500004	Right	Large	20.5-23.5"
B-237600004	Left	Large	20.5-23.5"
B-237500005	Right	XLarge	23.5-26.5"
B-237600005	Left	XLarge	23.5-26.5"

¹ Circumference: 6" above mid-patella



Medial View
showing Dynamic
Force Strap[®]