



Form Fit Tennis Elbow by Ossur

Description

The Form Fit Tennis Elbow has a patented ergonomically-shaped plastic design that eliminates the tourniquet effect associated with traditional tennis elbow brace designs.

The twin blue cushions provide a counterforce on were positioned on the inflamed tendon , reducing traction on the tendons and the tendon's attachment site where inflammation occurs.

The ribbed surfaces reduce vibration to the tendon, as well as increasing the comfort of the brace.

The Form Fit Tennis Elbow features a durable design and materials enable everyday use for both work and athletic activities.

Indications

- Tennis Elbow
- Golfers Elbow

Contact Us

Please feel free to contact for more information or to book an appointment

See more at: www.djoglobal.com/products/donjoy/quick-fit-wrist





CONRAD
BARNARD
ORTHOTICS & PROSTHETICS

E INFO@CONRADBARNARD.CO.ZA
T 033 342 4330 / 9464
C 083 301 4421
W CONRADBARNARD.CO.ZA
A 317 BULWER STREET, PIETERMARITZBURG, 3201